



A

LC/GC Mass Spectrometry



- Distinct ions possess a characteristic m/z and retention time.
- At each time, ions, with possibly overlapping m/z will be detected, and the ion count of ions with each m/z is recorded.
- Secondary MS (MS2) may be used to aid in ion identification or quantification.

Machine Output



Ion Identification

MS2

- m/z of daughter ions and parent ions are matched to possible peptides/metabolites.

Standards


- Pure standard with measured m/z and RT matched to current ions.
- Allows absolute quantification.

Peptide/metabolite-level analysis


- Based on the experimental design, estimate the main effects for each specie.
- The variance of species varies greatly.
- Individual measurements may be more variable than expected due to low signal quality.

B


Peptide-to-protein matching



Remove unsupported proteins and divergent peptides



Aggregate peptides to estimate protein abundance and variance



Downstream Applications